

## PLS TEST UPDATES EFFECTIVE IMMEDIATELY

---

### #907 Vitamin B6 (Pyridoxal 5-Phosphate)

Methodology: High Performance Liquid Chromatography  
Specimen Req: 1.0 mL (0.5 mL) Serum Non-SST or Plasma Heparinized (Lithium or Sodium)  
\*\*Protect from Light – Critical Frozen\*\*  
Separate samples must be submitted when multiple tests are ordered.  
Storage: Refrigerated 4 hours. Frozen 2 months.  
Set Up Days: Performed Sun-Sat. Reported 2-5 days.  
Performed At: ARUP #0080111

### #8140 MTHFR Mutation Detection PCR

**(PREVIOUS TEST NUMBER 9201: NOW DEACTIVATED)**

Methodology: Polymerase Chain Reaction  
Specimen Req: 3.0 mL Whole Blood EDTA.  
Storage: Refrigerated 1 week.  
Set Up Days: Thursdays. Reported 2-7 days.  
Performed At: PLS Molecular Department  
CPT: 83891, 83900, 83896 x 4, 83912  
Client Price: \$125.00

**PANEL COMPONENTS INCLUDE:**

**8141 MTHFR677T**  
**8142 MTHFRA1298C**  
**8139 INTERPRETATION**

### #1565 Troponin I

Methodology: Tosoh  
Specimen Req: 1.0 mL (0.5 mL) Plasma Heparinized (Serum No Longer Accepted).

### #2680 Exposure Panel

Methodology: Siemens Centaur Chemiluminescence  
Specimen Req: 3.0 mL Serum.  
Storage: Refrigerated 1 week.  
Set Up Days: Performed Mon-Fri.  
Performed At: PLS Omaha Chemistry  
CPT: 86703, 86803, 86706, 87340

**PANEL COMPONENTS:**

**680 HIV ½ Antibody**  
**945 Hepatitis C Antibody**  
**568 Hepatitis B Surface Antibody**  
**569 Hepatitis B Surface Antigen**

#457 Leukocyte Alkaline Phosphatase Stain (LAP Score)

**(PREVIOUS TEST NUMBER 223: NOW DEACTIVATED)**

Specimen Req: 5 mL (3 mL) Whole Blood Heparin.

Storage: Refrigerated 2 days.

**\*\*MUST BE RECEIVED BY PLS WITHIN 40 HOURS OF DRAW\*\***

**\*\*DO NOT SUBMIT FRIDAY – SUNDAY\*\***

Set Up Days: Performed Mon-Fri Only.

Performed At: UNMC Regional Lab

CPT: 85540

Client Price: \$140.00

#7695 Hantavirus Antibodies, IgG & IgM w/ Reflex to Confirmation

CPT: 86790 x 2; if reflexed add 86790 for each procedure.

Client Price: \$65.00; add \$60.00 for each procedure reflexed

**PANEL COMPONENTS MODIFIED:**

**7693 Hantavirus IgG**

**7694 Hantavirus IgM**

**8145 Sin Nombre IgG**

**8146 Sin Nombre IgM**

#1438 Respiratory Panel II

**PANEL COMPONENTS MODIFIED:**

**1449 Grass: Timothy**

**1462 Mold: A Alternata**

**1459 Peren: Cat Dander**

**1461 Peren: Cockroach**

**1460 Peren: Dog Dander**

**1457 Peren: Dust Mite DF**

**1458 Peren: Dust Mite DP**

**1473 Tree: Elm**

**1479 Weed: Ragweed**

#405 CK Isoenzymes

**PANEL COMPONENTS MODIFIED:**

**2130 Creatinine Kinase Total**

**1844 CK-MM**

**450 CK-MB**

**1845 CK-BB**

**8147 CK-Macro Type I**

**8148 CK-Macro Type II**

#8158 Phospholipid Antibodies, IgG, IgM, & IgA

**PANEL COMPONENTS MODIFIED:**

- 8149 Phosphatidylethanolamine Antibody IgG**
- 8150 Phosphatidylethanolamine Antibody IgM**
- 8151 Phosphatidylethanolamine Antibody IgA**
- 8152 Phosphatidylglycerol Antibody IgG**
- 8153 Phosphatidylglycerol Antibody IgM**
- 8154 Phosphatidylglycerol Antibody IgA**
- 8155 Phosphatidylinositol Antibody IgG**
- 8156 Phosphatidylinositol Antibody IgM**
- 8157 Phosphatidylinositol Antibody IgA**
- 7906 Phosphatidylserine Antibodies IgG, IgM, & IgA**
- 7272 Cardiolipin Antibody IgG**
- 7273 Cardiolipin Antibody IgM**
- 7274 Cardiolipin Antibody IgA**

IF YOU HAVE ANY QUESTIONS, PLEASE CONTACT KACEY MORELAND OR CRYSTAL WAYMIRE AT  
402-731-4145 OR 1-800-642-1117